


TRECVID 2012: instance.search results

Run ID: CeaList
 Processing type: automatic
 System training type: X (not specified)
 Condition: N (No IACC.1 *_meta.xml used)
 Priority: 2


Across 21 test topics (9048-9068)

Total relevant shots: 1232
 Total relevant shots returned: 202


Mean(prec. @ total relevant shots): 0.141
 Mean(average precision): 0.113



| Interpolated recall | precision | Precision at n shots |
|---------------------|-----------|----------------------|
| 0.0 | 0.4908 | 5 0.3333 |
| 0.1 | 0.3846 | 10 0.2857 |
| 0.2 | 0.2134 | 15 0.2571 |
| 0.3 | 0.1024 | 20 0.2333 |
| 0.4 | 0.0913 | 30 0.2016 |
| 0.5 | 0.0893 | 100 0.0871 |
| 0.6 | 0.0243 | 200 0.0479 |
| 0.7 | 0.0169 | 500 0.0192 |
| 0.8 | 0.0000 | 1000 0.0096 |
| 0.9 | 0.0000 | |
| 1.0 | 0.0000 | |



Run score (dot) versus median (---) versus best (box) by topic



Elapsed search time by topic