

TRECVID 2012: Semantic indexing (paired concepts) results


Run ID: P_A_TokyoTechCanon5_brn_5
 Processing type: Automatic
 Class: P - paired concepts
 Training type: A (only IACC data)
 Priority: 5

PLEASE NOTE: All of the measures below are based on assessment of a
 2-tiered random sampling (1-200@100%, 201-2000@10%) of the
 full submission pools and use of sample_eval.pl to infer the measures

Across 10 test concepts


Inferred total true shots: 2008
 Inferred true shots returned: 871

Mean(inferred average precision): 0.044



| Inferred interpolated recall precision | |
|--|-------|
| 0.0 | 0.298 |
| 0.1 | 0.115 |
| 0.2 | 0.090 |
| 0.3 | 0.072 |
| 0.4 | 0.050 |
| 0.5 | 0.028 |
| 0.6 | 0.018 |
| 0.7 | 0.005 |
| 0.8 | 0.000 |
| 0.9 | 0.000 |
| 1.0 | 0.000 |

| Inferred precision at depth n shots | |
|-------------------------------------|-------|
| 10 | 0.190 |
| 100 | 0.085 |
| 1000 | 0.055 |
| 2000 | 0.044 |



Run score (dot) versus median (---) versus best (box) by concept