


TRECVID 2012: instance.search results

Run ID: CeaList
 Processing type: automatic
 System training type: X (not specified)
 Condition: N (No IACC.1 *_meta.xml used)
 Priority: 3


Across 21 test topics (9048-9068)

Total relevant shots: 1232
 Total relevant shots returned: 166


Mean(prec. @ total relevant shots): 0.039
 Mean(average precision): 0.027



| Interpolated recall precision | | Precision at n shots | |
|-------------------------------|--------|----------------------|--------|
| 0.0 | 0.2420 | 5 | 0.1143 |
| 0.1 | 0.0675 | 10 | 0.0857 |
| 0.2 | 0.0491 | 15 | 0.0857 |
| 0.3 | 0.0322 | 20 | 0.0786 |
| 0.4 | 0.0140 | 30 | 0.0587 |
| 0.5 | 0.0040 | 100 | 0.0243 |
| 0.6 | 0.0000 | 200 | 0.0183 |
| 0.7 | 0.0000 | 500 | 0.0119 |
| 0.8 | 0.0000 | 1000 | 0.0079 |
| 0.9 | 0.0000 | | |
| 1.0 | 0.0000 | | |



Run score (dot) versus median (---) versus best (box) by topic



Elapsed search time by topic