

TRECVID 2013: Semantic indexing (paired concepts) results


Run ID: 13_P_A_PicSOM_P_5_5
 Processing type: Automatic
 Class: P - paired concepts
 Training type: A (only IACC data)
 Priority: 5

PLEASE NOTE: All of the measures below are based on assessment of a 2-tiered random sampling (1-200@100%, 201-2000@6.6%) of the full submission pools and use of sample_eval.pl to infer the measures

Across 10 test concepts

Inferred total true shots: 558
 Inferred true shots returned: 313

Mean(inferred average precision): 0.112



Inferred interpolated	
recall	precision
0.0	0.312
0.1	0.227
0.2	0.171
0.3	0.150
0.4	0.136
0.5	0.123
0.6	0.100
0.7	0.060
0.8	0.034
0.9	0.017
1.0	0.000

Inferred precision at depth n shots	
n	precision
10	0.250
100	0.123
1000	0.027
2000	0.016

